


MUTARS<sup>®</sup>-Münster



implantcast



**Distal Humerus**  
surgical technique



# MUTARS<sup>®</sup>-Münster

---

## Distal Humerus surgical technique

MUTARS<sup>®</sup> was developed in co-operation with  
Prof. Dr. W. Winkelmann,  
and Prof. Dr. G. Gosheger,  
Clinic and Polyclinic for General Orthopedics  
at the University Hospital of Münster, Germany.  
MUTARS<sup>®</sup> has been in successful clinical use since 1992

### Table of Contents


System overview.....	1
Assembling options/reconstruction length.....	2
Surgical technique .....	3
Implants .....	13
Instruments .....	15
Important medical information .....	19

**Nota Bene:** The described surgical technique is the suggested treatment for the uncomplicated procedure.  
In the final analysis the preferred treatment is that which addresses the needs of the individual patient.

**Copyright note:** MUTARS<sup>®</sup>, implavit<sup>®</sup> and implatan<sup>®</sup> are registered trademarks of implantcast GmbH. The use and copy of the content of this brochure are only allowed with prior permit given by the implantcast GmbH.



# System Overview






# MUTARS<sup>®</sup> Distal Humerus

assembling options  
(length in mm)

reconstruction length	humerus components				
	distal humerus	extension piece	connecting part 80 mm	extension piece	humerus screw
60	50				15
80	50	20			35
100	50	40			55
120	50	60			75
140	50		80		15 + 15
160	50	20	80		35 + 15
180	50	40	80		55 + 15
200	50	60	80		75 + 15
220	50	60	80	20	75 + 35
240	50	60	80	40	75 + 55
260	50	60	80	20 + 40	75 + 75

**Note:** Please notice that the amount of implants and instruments send with an individual shipment may differ from the information in the catalogue information of this brochure. Please make sure, during the preoperatively planning, that all necessary implants and instruments are available for the surgery.



**figure 1**

## Tumor resection

Measure the size of the resected amount of bone. Compare the length to the preoperative planning.

### Remark

In the case that the radius head is free of tumor it needn't be resected.

The minimum bone resection is 60mm.

## Preparation of the proximal ulna

Open the medullary cavity.  
Prepare the cavity with a sharp spoon (fig. 1) and an air-drill (fig. 2).

To implant the ulna anchorage in the exact depth, the entry has to be extended ventral at the Proc. Coronoideus and dorsal.

Rasp the medullary cavity with the MUTARS<sup>®</sup> rasp for ulna anchorage (fig. 3a and 3b).

### Remark

Risk of the via falsa with cortical perforation. Control in 2 plains under x-ray is recommended.



**figure 2**



**figure 3a**



**figure 3b**

## Implantation of the ulna anchorage

If a cemented implantation is planned, clean the ulnar cavity (fig. 4) and insert some cement.

Insert and impact the ulna anchorage cemented or cementless. Use either the straight or the cranked setting instrument for impaction (fig. 5).

Make again a concluding x-ray control in 2 plains.




figure 4




figure 5



## Screw fixation of the ulna anchorage

It is recommended to enhance the fixation of the ulna anchorage by adding a bone screw, both for the cementless and the **cemented** implantation.

Please follow the steps shown on the left. Drill with the 2mm drill (fig. 6a), measure the length of the 4 mm screw (fig. 6b) and insert the screw (fig. 6c).

### Remark

The use of a cancellous screw is preferable, because the use of a cortical screw can lead to a skin impingement.




figure 6a




figure 6b




figure 6c



## Humeral bone preparation

### Cementless use

Drill the medullary cavity with a humerus drill 1 mm smaller than the size of the preoperatively chosen humerus stem (fig. 7).

### Cemented use

Drill the medullary cavity with a humerus drill 2 mm larger than the size of the preoperatively chosen humerus stem (fig. 7).



figure 7

Prepare the medullary cavity with the medullary cavity reamer (fig. 8).



figure 8




figure 9

### Rasping of the humeral cavity

Assemble the humeral rasp of the appropriated size (see tables below), the extractor device, the humerus impactor and the sleeve. Lock the rasp on the humerus impactor by using the counter wrench.

### Remark

The use of a humeral rasp for a **cemented stem** is optional.

### Use of cementless stems

Use the humeral rasp (fig. 9) of the same size as the preoperatively chosen humerus stem (table 1).

Stem size	Rasp size
9mm	9mm
10mm	10mm
11mm	11mm
12mm	12mm
13mm	13mm

table 1

### Optional technique for the use of cemented stems

If you want to prepare for a cemented stem with the humeral rasp, please use the rasp which is 2 mm larger than the preoperatively chosen cemented humerus stem (fig. 9).

That will provide a cement mantle of 1mm thickness (table 2).

Stem size	Rasp size
8mm	10mm
9mm	11mm
10mm	12mm

table 2

Rasp the medullary cavity with the correct humeral rasp (fig. 10a and 10b). A carefully use of the mallet is recommended.

**Remark**

It is recommended to clean the rasp from bone chips during the rasping. Leave the humeral rasp in the bone for the trialing.

**Trial reduction**

Mount the distal humerus and the possibly used extension pieces (possible enlargement from 20 to 200mm; see table page2) onto the top of the rasp. Assemble the articulating mechanism by inserting the trial axle (fig. 11a).




figure 10a




figure 10b




figure 11a

figure 11b

## Implantation of the humeral stem

Impact the MUTARS<sup>®</sup> humerus stem (fig. 12).

Insert the stem of the same size as the rasp if a **cementless stem** is used. It is possible to protect the humerus stem against rotation by using a 4 mm cortical screw.

If a cemented implantation is planned insert the cement and use the **cemented stem** which is 2 mm smaller than the previously used drill or rasp.

Remove all instruments, especially during the cement hardening to prevent bending moments.




figure 12

## Implantation of the distal components

Combine the distal humerus on the humerus stem (fig. 13). If necessary extend with humerus extension pieces. Adjust the correct rotation position.

Lock the components with the corresponding humerus screw of the correct length (see table on page 2). Secure the components with the MUTARS<sup>®</sup> socket wrench small (fig. 14a).

Secure the assembly by using the counter instrument (fig. 14c).

Lock the humerus safety screw in the same way (fig. 14b).



figure 13



figure 14c



figure 14a



figure 14b



figure 15

Insert the ulna stop with the setting instrument for ulna anchorage (fig. 15).

The ulna stop must entirely close the hole of the distal humerus to ensure a free run of the joint.



figure 16

### Final trial reduction

Connect the distal humerus to the ulnar anchorage by inserting the articulating axle (fig. 16).

## Locking of the articulating axle

To cover the articulating mechanism and to protect the axle on both sides the locking screws are inserted (fig. 17a).

Therefore the MUTARS<sup>®</sup> socket wrench small is used (fig. 17b and 17c).



figure 17a




figure 17b



figure 17c



# MUTARS<sup>®</sup> Distal Humerus

## IMPLANTS

### MUTARS<sup>®</sup> Distal humerus 50 mm incl. axle, safety and 2 lock screws humerus cap

mat.: implatan<sup>®</sup>; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3  
 axle CoCrMo forged alloy according to  
 DIN ISO 5832/12  
 5250-0000



### MUTARS<sup>®</sup> ulna anchorage cementless

mat.: implatan<sup>®</sup>; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3  
 with cpTi and HA-coating  
 5250-1015 left  
 5250-1020 right



### MUTARS<sup>®</sup> ulna anchorage cemented

mat.: implavit<sup>®</sup>; CoCrMo-casting alloy according to  
 DIN ISO 5832/4  
 5250-1005 left  
 5250-1010 right



### MUTARS<sup>®</sup> ulna stop

mat.: UHMWPE according to DIN ISO 5834/2  
 5250-1100

### locking bone screw 4mm

mat.: implatan<sup>®</sup>; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3  
 5793-4026 26mm  
 5793-4028 28mm  
 5793-4030 30mm  
 5793-4032 32mm  
 5793-4034 34mm



### MUTARS<sup>®</sup> attachment tube

mat.: polyethylenterephthalat  
 5900-0300 35 mm  
 5900-0310 55 mm





### IMPLANTS

#### MUTARS® humerus screw

mat.: *implatan*®;  $TiAl_6V_4$  according to DIN  
ISO 5832/3

5230-0015	M8x15 mm
5230-0035	M8x35 mm
5230-0055	M8x55 mm
5230-0075	M8x75 mm



#### MUTARS® humerus stem cementless

mat.: *implatan*®;  $TiAl_6V_4$  according to DIN  
ISO 5832/3

5240-0709	9 mm
5240-0710	10 mm
5240-0711	11 mm
5240-0712	12 mm
5240-0713	13 mm



#### MUTARS® humerus stem cemented

mat.: *implavit*®; CoCrMo-casting alloy  
according to DIN ISO 5832/4

5240-0408	8 mm
5240-0409	9 mm
5240-0410	10 mm
5240-0411	11 mm
5240-0412	12 mm

**special stem sizes are available on request.**



#### MUTARS® humerus extension piece

mat.: *implatan*®;  $TiAl_6V_4$  according to DIN  
ISO 5832/3

5220-0020	20 mm
5220-0040	40 mm
5220-0060	60 mm






# MUTARS® Distal Humerus

---

## INSTRUMENTE



**MUTARS® humerus instrument tray**  
7999-5200



**MUTARS® Distal humerus instrument tray**  
7999-5150



**MUTARS Humerus trial component tray**  
7999-5202



**INSTRUMENTS**

**MUTARS® extractor device**  
7220-0000



**MUTARS® socket wrench small**  
7608-1010



**MUTARS® medullary cavity reamer**  
7760-0500



**MUTARS® humerus drill**  
7 mm      7630-0007  
8 mm      7630-0008  
9 mm      7630-0009  
10 mm     7630-0010  
11 mm     7630-0011  
12 mm     7630-0012



**MUTARS® humeral rasp**  
9 mm      7770-0709  
10 mm     7770-0710  
11 mm     7770-0711  
12 mm     7770-0712  
13 mm     7770-0713



**MUTARS® humerus impactor**  
7710-0000



## INSTRUMENTS



**MUTARS® humerus  
impact and extract sleeve**  
7721-0000



**MUTARS®  
wrench for humerus**  
7420-0001



**MUTARS® humerus  
wrench for cap/  
counter instrument**  
7710-0001



**MUTARS® humerus trial cap**  
7710-1000 small  
7710-1005 medium  
7710-1010 large



**MUTARS® rasp for ulna anchorage**  
left 7420-0006  
right 7420-0007



**MUTARS® instrument for  
ulna anchorage straight**  
7420-0013



**MUTARS® instrument for  
ulna anchorage cranked**  
7420-0014



**MUTARS® trial axle  
for distal humerus**  
7420-0015



**depth measuring instrument small**  
0270-1015



**drill 2 mm**  
7520-0000



**hexagon head screw driver 2,5 mm**  
7608-1001



**ball reamer 4 mm**  
7700-2604



## INSTRUMENTS

### MUTARS® humerus trial cap with taper

- 7710-1200 small
- 7710-1205 medium
- 7710-1210 large



### MUTARS® humerus trial head

7710-1252



### MUTARS® humerus trial extension piece

- 7710-0020 20 mm
- 7710-0040 40 mm
- 7710-0060 60 mm



### MUTARS® humerus trial reducer

- 7710-2100 10 mm
- 7710-2101 100 mm



### MUTARS® humerus trial connecting part 80 mm

7710-2180



### MUTARS® humerus trial screw

- 7710-2315 M8x15 mm
- 7710-2335 M8x35 mm
- 7710-2355 M8x55 mm
- 7710-2375 M8x75 mm





# IMPORTANT MEDICAL INFORMATION

## Note

The MUTARS® system is a successful therapy, freeing e.g. tumor patients from pain and restricted mobility.

The system's main objectives are pain reduction and the restoration of physiological functions. Suitable patients should be selected according to the following conditions:

- 1) Patients disposing of adequate bone quality and with a bone structure that is likely to be preserved.
- 2) Patients, whose anatomic features allow for an adequate implant size for the prospective loading and degree of activity.
- 3) Patients who are willing and able to follow their physician's directions, especially with respect to the necessary total or partly stress reduction on the implant during the postoperative period.

The largest possible stem size is to be selected from the MUTARS® system (especially for obese patients). Patients must be warned of the consequences of excessive weight-bearing, sport participation and any activity causing excessive strain or impingement on the implanted prosthesis.

## Used materials

The MUTARS® implants consist of a cast CoCrMo-alloy (ISO 5832/4) or the titanium alloy TiAl6V4 (ISO 5832/3). The PE-components contain UHMW-PE (ISO 5834/2). Most of the instruments are made of acid-resistant stainless steel.

## Indications

Adequate patient selection as well as a profound surgical analysis of the case are the basis of the whole surgical procedure. Careful preoperative planning and a precise surgical technique are necessary to obtain optimal results. In order to minimize the danger of postoperative complications different factors must be considered, i.e. the anatomical stress situation, the soft tissue basis and the planned component alignment. Generally, a prosthesis is only to be implanted in patients with fully-grown skeletons.

To restore the anatomical function of the skeleton it may be necessary to repair and/or support/stabilize a traumatised or otherwise affected bone segment, to fuse it with other fragments or replace it by a prosthesis.

The treatment of fractures, pseudarthrosis, arthrosis and similar diseases with the MUTARS® system can either be performed as an initial or a follow-up surgery, each with its corresponding surgical technique.

The MUTARS® system is mainly implanted in cases of major bone defects, e.g. after bone tumor resections.

The use of a modular prosthesis is often the consequence of a tumor resection. In the case of primary tumors it is necessary to perform an

extensive resection (Enneking-method), extending into healthy tissue, in order to provide for an adequate surgical treatment of the disease. If this is not possible, other steps must be discussed, such as e.g. an amputation. The tumor system is not intended to support an intraligamentary/marginal and - with respect to the stage of the disease - in such cases unsuitable therapy.

In the case of bone metastases the indication depends on the entire patient status. Whenever a skeletal segment can take no more strain and stabilizing fusions steps of an osteosynthesis are not possible, the tumor implant system can restore the function of the segment within a short period of time. This increases the patient's life quality considerably. Nevertheless, this indication must be reconsidered in the case of a multiple invasion of the bone in which a remobilization of the patient cannot be expected.

## Contraindications

The main contraindications are bacterial infections, soft tissue defects due to irradiation, expected bone growth as well as - under certain conditions - better alternatives as e.g. a resection arthrodesis for infants and growing young persons or a temporary prosthesis. The major contraindication is when dimension and localization of the tumor make an extensive resection impossible. Other contraindications include:

- 1) Anatomical conditions which do not allow for or will not maintain sufficient bony support of the implant or which do not allow for an adequate implant size.

In general i.e.:

- a) Insufficient blood supply caused by preceding surgeries or vessels affected by alcohol abuse etc.,
  - b) insufficient quantity and quality of bone material due to osteoporosis, adipositas etc.,
  - c) infections or other causes leading to reduced stability of the implant fixation.
- 2) Any mental or neurological conditions affecting the patient's will to follow restrictions in activity, especially during the post-operative healing process. These could be drug abuse, mental illness, senility and general neurological limitations.
  - 3) Conditions leading to extreme stress on the implants, such as myopathies, multiple arthropathies etc.

Contraindications can be absolute or relative and must be carefully considered with respect to the whole patient status and the prognoses of possible alternative therapies such as e.g. a conservative treatment, an arthrodesis etc.

## Possible adverse effects:

- 1) Device component loosening, distortion or fracture. Normally these effects are caused by one or several of the mentioned factors, listed above and below under contraindications and warnings.
- 2) Migration, subluxation and rotation of the implant, flexion contraction, reduced mobility, increased or decreased leg length, component loosening or bone wear and ligamentary laxity.
- 3) Fractures of the tibia, femur, patella and humerus.
- 4) Acute postoperative wound infection, severe sepsis and/or low-grade synovitis.
- 5) Neuropathy
- 6) Cardiovascular disorders: wound-haematoma, thrombosis and embolism (e.g. venous thrombosis and pulmonary embolism)
- 7) Tissue reactions: phagocytal reactions, foreign body reactions and myositis ossificans. These reactions especially apply to male patients with hypotrophic arthrosis, preoperatively limited mobility and/or preceded myositis. The risk of a myositis ossificans is increased by preceding surgeries or acute infections.
- 8) Trochanteric pseudarthrosis: generally related to early stress and/or insufficient fixation in the case of a transtrochanteric surgical path.

## Warning and Precautions

Implant loosening, bending, fissure and/or breakage and other complications can occur if the following instructions and warnings are not considered and followed

### Preoperative:

- 1) In every surgery a sufficiently wide range of implant sizes must be present. The decision, whether cementation is to be performed or not, must be taken in advance. The preoperatively chosen implant, as well as the next bigger and smaller sizes must be prepared. Before insertion, the implant must be carefully checked to make sure that there is no damage or modification and that the correct size has been selected.
- 2) The implants are to be handled with care at all times, in order to avoid damage of the prepared implant surface. Cutting, bending or scratching of the component surfaces can considerably reduce their stability and resistancy against fatigue and wear. Even not directly visible defects can cause stress conditions within the implant, which can - because of the dynamic stress within the body - possibly lead to implant failure. If the preoperative observation shows that the modularity of the system can not fit the patient the use of a customised implant is necessary.
- 3) Allergies and other reactions to implanted materials should be considered, tested (if indicated) and excluded preoperatively, even if very uncommon.
- 4) The introducing instruments must correspond to the implant and must therefore belong to the MUTARS® system.
- 5) A description of the surgical technique is available from the manufacturer. In order to obtain best possible results the surgeon must be familiar with the recommended surgical techniques for this system and its proper use.

### Intraoperative

- 1) Adequate and durable component support, obtained through cementation and/or bone material, as well as the correct selection of the component size are of vital importance.
- 2) Whenever a stem cementation is performed the entire stem must be cemented right up to the stem plate. During the process of cement hardening any repositioning of the implant components should be avoided.
- 3) After insertion of the stem, its plate must be flush against the resected bone. It is important to resect the bone plane, horizontal to the medullary canal.
- 4) For cementless tibia and femoral stems the use of our special MUTARS® rasps is mandatory.
- 5) Correct axial and rotational alignment of the implant is of major importance. Otherwise subluxation, dislocation and/or breakage of the prosthesis may occur. Special attention should be directed to cases with curved stems, since fixation might be achieved unplanned by a rotation of the implant during the insertion of the stem. In this case the implant-bone interface is insufficient.
- 6) In cases of congenital dysplastic coxarthropathy special care must be directed to the avoidance of a sciatic nerve paralysis. Moreover the fact must be considered that the medullary canal is often extremely narrow and straight, so that extremely small, straight femoral prostheses are necessary. Nevertheless the standard size should be applied whenever possible. Please consider that in these cases the original acetabulum is formed only rudimentarily and very narrow. Because of its anatomical biomechanical unreliability the acetabulum should not be used as implant bearing for the acetabulum component of the prosthesis.

7) From the technical point of view the performance of a revision surgery after preceding primary surgery is extremely demanding and critical. Common mistakes are: wrong surgical access, insufficient bone identification and mobilization, insufficient removal of ectophytic bone material or unprecise component positioning. Postoperative instability as well as e

xtrême blood loss can be the consequences. Altogether longer operating times, increased blood loss and the risk of pulmonary embolism and wound haematoma must be taken into consideration in cases of revision surgery.

8) Conus surfaces must be thoroughly cleaned and dried before attaching the fitting component. Any unremoved particle can cause extreme friction and wear.

9) Implants whose conus has been attached to an endohead before should not be reused, since the conus interface has adapted itself to the former endohead. A new endohead would therefore not fit properly.

10) After bar screw tightening with the MUTARS® swing wrench and articulated MUTARS® engineers' wrench SW 24, the bar screw should be countered in order to obtain the necessary fixation.

11) To avoid damage of the threads the bar screws should always be tightened completely.

## Postoperative

1) Postoperative patient care, detailed instructions and warnings by the physician are of major importance.

To enhance the healing process an external support of the operated leg is recommended for a limited period of time.

2) Active and passive movement must be carried out with extreme caution.

3) Postoperative therapy should support the healing process and prevent the leg from being submitted to excessive stress.

4) Patients are to be reminded repeatedly of the necessity to reduce activity as recommended by their physician.

## Special user information

Never reuse implants which have already been implanted or removed, even if they appear undamaged (danger of implant breakage due to internal material fatigue).

## Packaging and labelling

Each of the MUTARS® implants and instruments is packaged separately. The packaging of the not steril products is suitable for steam and Ethyloxyl-sterilization. The MUTARS® PE-components are supplied sterile. They should only be accepted by hospitals and physicians if supplied in their original packaging and with an original label.

For reasons of safety, protection and identification all implants should always be kept in cool and dry environment in their unopened packaging.

## Sterilization

The implants of the MUTARS® implant system have been sterilized by gamma-radiation (min. 25kGy) and are supplied in protective covers. Please always check the packages for perforation or other damage prior to surgery.

Resterilization of PE-components is not permitted.

For further information please refer to:



implantcast GmbH

Lüneburger Schanze 26

D-21614 Buxtehude, Germany

phone : +49 4161 744 0

fax : +49 4161 744 200

e-mail: [info@implantcast.de](mailto:info@implantcast.de)

internet : [www.implantcast.de](http://www.implantcast.de)



Technical modifications are subject to change without notice. All implants have been developed, manufactured and tested according to latest technical standards. No part of this document may be duplicated without prior consent of implantcast



implantcast GmbH  
Lüneburger Schanze 26  
D-21614 Buxtehude  
Germany  
phone: +49 4161 744-0  
fax: +49 4161 744-200  
e-mail: [info@implantcast.de](mailto:info@implantcast.de)  0123  
internet: [www.implantcast.de](http://www.implantcast.de)

Your local distributor:

